



Holder Type	Description	Applications	Maximum Size mm (in) WxHxD	Exposure Size mm (in) WxH	Instrument	Rack Capacity
 SL-3T Single exposure window with spring clip back Part Number 191163900	General purpose holder for exposure of thin specimens	Textiles, plastic film, automotive interior	69 x 145 x 3 (2.7 x 5.7 x 0.1)	50 x 121 (2.0 x 4.7)	Ci3000+ Ci4000 Ci5000	20 68 111
 SL-3T with Glass Single exposure window with glass and adjustable back Part Number 07303900	Exposure through secondary filter to closely simulate end-use environment	Textiles, paper, plastic film, carpet, automotive interior	69 x 145 x 15 (2.7 x 5.7 x 0.6)	50 x 121 (2.0 x 4.7)	Ci3000+ Ci4000 Ci5000	20 68 111
 CD-3T Three exposure windows with spring clip back Part Number 20215700	General purpose holder with three exposure windows for thin specimens	Textiles, paper, plastic film, automotive interior	69 x 145 x 3 (2.7 x 5.7 x 0.1)	3 windows: 38 x 50 (1.5 x 2.0)	Ci3000+ Ci4000 Ci5000	20 68 111
 CD-3T with Glass Three exposure windows with glass, spring clip back Part Number 07303800	Exposure through secondary filter to closely simulate end-use environment	Textiles, paper, plastic film, wood, automotive interior	69 x 145 x 15 (2.7 x 5.7 x 0.6)	3 windows: 38 x 50 (1.5 x 2.0)	Ci3000+ Ci4000 Ci5000	20 68 111
 CD-2W Thick carpet, hinged with support Part Number 07255500	Non-compressing exposure of thick, non-rigid materials	Carpet, foam, foam-backed materials	71 x 145 x 12 (2.8 x 5.7 x 0.5)	60 x 66 (2.3 x 2.6)	Ci3000+ Ci4000 Ci5000	20 68 111
 WPTC-3T Part Number 06150400	Single window opening to expose wallpaper or thick carpeting	Carpet, foam, foam-backed materials, patterned materials	165 x 146 x 12 (6.5 x 5.7 x 0.5)	131 x 100 (5.1 x 3.9)	Ci3000+ Ci4000 Ci5000	8 27 39
 TEX-3T with Mask Single exposure window with mask, adjustable Part Number 19186700	Holder with masked exposure and adjustable back to allow for thick sample	Textiles, foam, foam-backed materials	45 x 134 x 12 (0.7 x 5.3 x 0.5)	19 x 119 (.07 x 4.7)	Ci3000+ Ci4000 Ci5000	29 104 170
 Reference Chip Polystyrene reference chip holder Part Number 19183400	Holds a single chip of polystyrene automotive reference material	Polystyrene reference chips	50 x 88 x 2 (2.0 x 3.5 x 0.1)	43 x 82 (2.8 x 3.6)	Ci3000+ Ci4000 Ci5000	16 54 93

This chart is a representative sample of specimen holders available for a variety of Atlas instruments. For specific information about specimen holders that best meet your needs, please contact your local Atlas representative.



Specimen Holders

Rack Capacity	Instrument	Exposure Size mm (in) WxH	Maximum Size mm (in) WxHxD	Applications	Description
12 41 71	Ci3000+ Ci4000 Ci5000	101 x 146 (4.0 x 5.7)	104 x 155 x 12 (4.1 x 6.1 x 0.5)	Coatings, rigid plastic, wood	Holds 4 x 6 panels for exposure tests of coatings
16 56 87	Ci3000+ Ci4000 Ci5000	76 x 146 (3.0 x 5.7)	76 x 152 x 9 (3.0 x 6.0 x 0.37)	Coatings, rigid plastic, wood	Holds 3 x 6 panels for exposure tests of coatings
9 35 51	Ci3000+ Ci4000 Ci5000	119 x 119 (4.7 x 4.7)	127 x 138 x 9 (5.0 x 5.4 x 0.3)	Rigid plastic, roofing material, solar panels, wood	Holds large, self supporting material for exposure
20 65 111	Ci3000+ Ci4000 Ci5000	50 x 121 (2.0 x 4.7)	69 x 101 x 43 (2.7 x 4.0 x 1.7)	Bottles, labels, printing inks, adhesives, liquids, pills	Holds small bottles or test tubes for exposure of packaging and pharmaceuticals
15 50 84	Ci3000+ Ci4000 Ci5000	76 x 125 (3.0 x 4.9)	77 x 144 x 3 (3.0 x 5.6 x 0.1)	Plastics	Holds vertical specimen bars
15 50 75	Ci3000+ Ci4000 Ci5000	71 x 121 (2.8 x 4.7)	85 x 145 x 3 (3.7 x 5.7 x 0.15)	Plastics	Holds an 85 mm tensile bar
20 68 111	Ci3000+ Ci4000 Ci5000	56 x 127 (2.2 x 5.0)	55 x 137 x 5 (2.2 x 5.4 x 0.2)	Plastics	Holds specimens of varying sizes and shapes: bars, discs
21 68 111	Ci3000+ Ci4000 Ci5000	39 x 138 (1.5 x 5.4)	50 x 151 x 3 (2.0 x 5.9 x 0.1)	35 mm slides, rigid discs, plaques	Holds 50 x 50 mm slides for exposure
14 41 71	Ci3000+ Ci4000 Ci5000	101 x 92 (4.0 x 3.6)	101 x 101 x 10 (4 x 4 x 0.4)	Automotive or building glass	Holds varying widths of glass attaching at top and bottom for exposure

Holder Type

Panel Holder
Part Number 19210200

Panel Holder
Part Number 19188501

Solar Panel Holder
Part Number 19190400

Adjustable Bottle Holder
Part Number 19178100

Drop-in Specimen Bar Holder
Part Number 19184600

Tensile Bar Holder with Spring Clip Back
Part Number 19212100

Adjustable Specimen Holder
Part Number 19210600

Slide Holder
Part Number 19195800

Glass Holder
Part Number 19181900





Holder Type	Description	Applications	Maximum Size mm (in)	Exposure Size mm (in)	Instrument	Rack Capacity
 Regular Specimen Holder Part Number 56075142	Holds specimen up to 3 mm thick	Textiles, plastics, coatings, paper	135 x 45 (5.3 x 1.8)	121 x 35 (4.8 x 1.4)	Xenotest® Alpha Xenotest® 150 S+	11
 Specimen Holder for Thick Specimens Part Number 56077987	Holds specimen up to 15 mm thick	Carpets, foam, foam-backed materials, thick panels	135 x 45 (5.3 x 1.8)	121 x 35 (4.8 x 1.4)	Xenotest® Alpha Xenotest® 150 S+	11
 Specimen Holder for Blue Scale Part Number 56050873	Holder for Blue Scale fabric during weathering tests	Blue Scale fabric	135 x 45 (5.3 x 1.8)	—	Xenotest® Alpha Xenotest® 150 S+	only 1 needed (replaces one holder)
 Standard Specimen Holder (including back plate) Part Number 56076543	Holds specimen up to 10 mm thick	Textiles, plastics, coatings, paper, window profiles	310 x 80 (12.2 x 3.2)	295 x 70 (11.6 x 2.8)	Xenotest® Beta+	16
 Special Specimen Holder, 3 Segments Part Number 56076699	Holds 3 specimens up to 10 mm thick, fixed by a back plate and a cover mask	Textiles, plastics, coatings, paper	100 x 68 (3.9 x 2.7)	90 x 58 (3.5 x 2.3)	Xenotest® Beta+	16
 Special Specimen Holder, 2 Segments Part Number 56076700	Holds 2 specimens up to 10 mm thick, fixed by 2 cover masks	Textiles, plastics, coatings, paper	135 x 45 (5.3 x 1.8)	121 x 35 (4.8 x 1.4)	Xenotest® Beta+	28
 Specimen Holder for Blue Scale (Special holder kit, 2 Segments) Part Number 56078036	Holds Blue Scale fabric for weathering tests	Blue Scale fabric	135 x 45 (5.3 x 1.8)	—	Xenotest® Beta+	only 1 needed (replaces one holder)
 Set of five specimen trays for 7 specimens Part Number 56079259	General purpose holder for various materials	Plastics, coatings, paper	135 x 45 (5.3 x 1.8)	121 x 35 (4.8 x 1.4)	SUNTEST XXL/XXL+	34
 Set of two specimen trays for four specimens Part Number 56079261	General purpose holder for various materials	Plastics, coatings, paper	310 x 80 (12.2 x 3.2)	295 x 70 (11.6 x 2.8)	SUNTEST XXL/XXL+	8

This chart is a representative sample of specimen holders available for a variety of Atlas instruments. For specific information about specimen holders that best meet your needs, please contact your local Atlas representative.